



GLP Compliance in Digital Pathology

Curtis Adams Ph.D.

Sr. Product Manager – Life Sciences

curtis@aperio.com

June 30, 2011

www.aperio.com

Agenda

- General Considerations
 - **General Laboratory Practices**
 - history
 - agencies
 - guidelines
 - observations
 - **Comprehensive Solution to Regulatory Guidelines**
 - instrument qualification
 - software validation
- Advantages for Compliance

Agenda

- General Considerations
 - **General Laboratory Practices**
 - history
 - agencies
 - guidelines
 - observations
 - Comprehensive Solution to Regulatory Guidelines
 - instrument qualification
 - software validation
- Advantages for Compliance

History of GLP

- **GLP Regulations** describe the minimal standards for conducting nonclinical laboratory studies the support applications for research or marketing permits for products regulated by FDA or EPA such as **human drugs**, food additives, medical devices, biological products, and pesticide products.
- **Origins:** GLP regulatory mandate was first promulgated in 1978 by the US-FDA and published in the Federal Register 43 FR 59985-60020. Since then the Organisation for Economic Co-operation and Development (OECD) helped promulgate it to many countries to place in their national regulations.



OECD Member Countries

- Australia
- Austria
- Belgium
- Canada
- Czech Republic
- Denmark
- Finland
- France
- Germany
- Greece
- Hungary
- Iceland
- Ireland
- Italy
- Japan
- S. Korea
- Luxembourg
- Mexico
- Netherlands
- New Zealand
- Norway
- Poland
- Portugal
- Slovak Republic
- Spain
- Sweden
- Switzerland
- Turkey
- United Kingdom
- United States

Alphabet soup of regulations



GLP (21 CFR § 58)

ERES (21 CFR § 11)

GCP (21 CFR § 312, etc)

GMP (21 CFR §210 and 211)

Compliance Policy Guides

Guidances for Industry



OECD Principles of Good Laboratory Practice

Guidance Documents

ERES (Annex 11, “Computerised Systems”)



USP <1058>

ISPE GAMP 4 and 5

HIPAA (45 CFR § 160 and 164)

Q9 (Quality Risk Management)

Q10 (Quality System)

E6 (GCP)

E7-E11 (Clinical Trials)



GLP (MHLW Ordinance No. 21)

GCP (Ordinances and Notifications)

ERES (PFSB Notification 0401022)

Overview of ERES and predicate GLP and GCP requirements

Stated or Implied Requirements from ERES regulations, GLP, and GCP Requirement	USA 21 CFR Part 11 “Electronic Records; Electronic Signatures”	EU Annex 11 “Computerized Systems”	Japan PFSB Notification, “Use of Electromagnetic Records and Electronic Signatures...”	Predicate GLP and GCP regulations (USA, EU, JP)
Risk management system must be in place		➡		
Instruments, software and systems must be validated	➡	➡	➡	➡
System components must be inventoried		➡		➡
Equipment must be adequately tested, calibrated, and standardized	➡	➡	➡	➡
System generates accurate and complete copies of records for inspection	➡		➡	
Records must be protected	➡	➡	➡	
System access must be limited	➡	➡	➡	
System must identify who created records	➡	➡	➡	➡
Audit trails must track who created/changed/deleted records	➡	➡	➡	➡
Authority checks must be in place	➡	➡	➡	
Device / terminal checks must be in place	➡			
Additional controls must be placed on open systems	➡			
Electronic signature manifestations indicate who / when / why	➡	➡	➡	
Electronic signature irrevocably linked to record	➡	➡	➡	
Electronic signature components and controls are enforced	➡	➡	➡	
System alerts of bad login and bad e-sig attempts	➡	➡		
Data backups are taking place routinely		➡	➡	➡

Why CFR Section 21?

- *Code Of Federal Regulations (CFR)*
- The final regulations published in the *Federal Register* (daily published record of proposed rules, final rules, meeting notices, etc.) are collected in the *CFR*. The *CFR* is divided into 50 titles which represent broad areas subject to Federal regulations. The FDA's portion of the *CFR* interprets the *Federal Food, Drug and Cosmetic Act* and related statutes. *Section 21 of the CFR* contains all regulations pertaining to food and drugs. The regulations document all actions of all drug sponsors that are required under Federal law.

Source:  U.S. Food and Drug Administration

CFR 21 part 58 – General Laboratory Practices

Sec. 58.1 Scope.

- (a) This part prescribes good laboratory practices for conducting nonclinical laboratory studies that support or are intended to support applications for research or marketing permits for products regulated by the Food and Drug Administration, including food and color additives, animal food additives, **human and animal drugs**, medical devices for human use, biological products, and electronic products. Compliance with this part is intended to assure the quality and integrity of the safety data filed pursuant to sections 406, 408, 409, 502, 503, 505, 506, 510, 512-516, 518-520, 721, and 801 of the Federal Food, Drug, and Cosmetic Act and sections 351 and 354-360F of the Public Health Service Act.
- (b) References in this part to regulatory sections of the Code of Federal Regulations are to chapter I of title 21, unless otherwise noted.

[43 FR 60013, Dec. 22, 1978, as amended at 52 FR 33779, Sept. 4, 1987; 64 FR 399, Jan. 5, 1999]

CRF 21 part 58.3(k) definitions

(k) **Raw data** means any laboratory worksheets, records, memoranda, notes, or exact copies thereof, that are the result of original observations and activities of a nonclinical laboratory study and are necessary for the reconstruction and evaluation of the report of that study. In the event that exact transcripts of raw data have been prepared (e.g., tapes which have been transcribed verbatim, dated, and verified accurate by signature), **the exact copy or exact transcript may be substituted for the original source as raw data.** *Raw data* may include photographs, microfilm or microfiche copies, computer printouts, magnetic media, including dictated observations, and recorded data from automated instruments.

Source:  U.S. Food and Drug Administration

Required SOPs

- **Sec. 58.81** Standard operating procedures. (a) A testing facility shall have standard operating procedures in writing setting forth nonclinical laboratory study methods that management is satisfied are adequate to insure the quality and integrity of the data generated in the course of a study. All deviations in a study from standard operating procedures shall be authorized by the study director and shall be documented in the raw data. Significant changes in established standard operating procedures shall be properly authorized in writing by management.
- (b) Standard operating procedures shall be established for, but not limited to, the following:
 - (1) Animal room preparation.
 - (2) Animal care.
 - (3) Receipt, identification, storage, handling, mixing, and method of sampling of the test and control articles.
 - (4) Test system observations.
 - (5) Laboratory tests.
 - (6) Handling of animals found moribund or dead during study.
 - (7) Necropsy of animals or postmortem examination of animals.
 - (8) Collection and identification of specimens.
 - (9) **Histopathology.**
 - (10) Data handling, storage, and retrieval.
 - (11) Maintenance and calibration of equipment.
 - (12) Transfer, proper placement, and identification of animals.
- (c) Each laboratory area shall have immediately available laboratory manuals and standard operating procedures relative to the laboratory procedures being performed. Published literature may be used as a supplement to standard operating procedures.
- (d) A historical file of standard operating procedures, and all revisions thereof, including the dates of such revisions, shall be maintained.
- [43 FR 60013, Dec. 22, 1978, as amended at 52 FR 33780, Sept. 4, 1987]

Source:  **U.S. Food and Drug Administration**

Comprehensive Approach : Raw data: digital slide or glass slide?

Interpretations of virtual images have been repeatedly shown to be as accurate as interpreting glass slides...

Gold Standard – glass slides

75%

Agreement with peers on
“clinical significance”

90%

Agreement with self on
“clinical significance”

Histopathology 2007, 50, 266-173

A randomized controlled trial of the diagnostic accuracy of internet based telepathology compared with conventional microscopy

P. Furness

“No significant difference in diagnostic accuracy could be detected between the diagnoses proffered on the basis of virtual slides and conventional slides

Histopathology 2002, 41, 91-109

Telepathology: current status and future prospects in diagnostic histopathology

S. S. Cross, T. Dennis & R. D. Start


“The review concludes that all the necessary technology for telepathology is available and there is strong published evidence for a diagnostic accuracy comparable with glass slide diagnosis...”

Aperio’s FDA clearances

... which makes the raw data (image files) and metadata (annotations and processing) subject to ERES controls

Comprehensive Approach: Raw data: digital slide or glass slide?

- Electronic revolution: ERES, EHR, ERSR, etc
- Logical solution to regulatory requirements



- Electronic records enable integration across databases with inherent security, tracking, archiving/retrieval and back-ups.

(Further detail in *Toxicologic Pathology* 2007, 450-455, D. Tuomari et al, "Society of Toxicologic Pathology Position Paper on Pathology Image Data: Compliance with 21 CFR Parts 58 and 11")

Society of Toxicologic Pathology conclusions

The Society of Toxicologic Pathology recommends that **images used for data generation** (e.g., the basis of a diagnosis or morphometric analysis) **are raw data**, and in contrast, that images not used for data generation are illustrative images that are not raw data. Based on current technologies and practices, any image used for data generation, becomes raw data at the time of data generation and at that time an image print or the electronic image record must be authenticated by specific annotation indicating when and who used that image for data generation, and also that image raw data must be archived.

Source: *Toxicologic Pathology* 2007, 450-455, D. Tuomari et al,
“Society of Toxicologic Pathology Position Paper on Pathology Image Data:
Compliance with 21 CFR Parts 58 and 11

Form FDA 483 – Inspectional Observations

- Disclaimer – Form FDA 438 contains the observations of the inspector and does not necessarily “represent a final Agency determination regarding your compliance”.
- FDA 483 observations should listed in order of significance and may include previous observations that have not been corrected.
- Only those observations directly linked to a violation of regulations are typically included. Suggestions, guidance or other comments are typically not included.
- FDA publishes select 483’s on their website:
<http://www.fda.gov/AboutFDA/CentersOffices/ORA/ORAElectronicReadingRoom/default.htm>

Consequences of not validating GLP lab systems

FDA inspection observations have included:

- “Software... has not been fully validated for its intended use according to an established protocol. Electronic records are used, but they do not meet requirements to ensure that they are trustworthy, reliable, and generally equivalent to paper records”
- “There was a failure to check for accuracy of the inputs to and outputs from the TotalChrom Data Acquisition System, which is used to run your firm’s HPLC instruments ...”




Agenda


- General Considerations
 - **General Laboratory Practices**
 - history
 - agencies
 - guidelines
 - observations
 - **Comprehensive Solution to Regulatory Guidelines**
 - instrument qualification
 - software validation
- Advantages for Compliance

Digital Pathology Components

Slide scanners to digitize whole slides at typical study volumes



Secure data repository for storage, retrieval, analysis, plotting, reporting, archiving, etc



On-site services for GLP validation, international multisite integration etc






Image analysis for whole slide, automated, objective data




Tissue Microarray Lab for biomarker discovery



Qualification vs Validation

- Qualification : the verification that **an instrument** is performing under predetermined specifications.
- Validation: the process of evaluating the performance of a specific measuring procedure and checking that the performance meets certain preset criteria.



Definitions: IQ, OQ, and PQ

	Definition	Examples
IQ	<p>INSTALLATION QUALIFICATION: <i>Documented verification that system components and support systems have been installed correctly and completely, in accordance with manufacturer and customer requirements, government regulations, and industry standards.</i></p>	<ul style="list-style-type: none"> • Correct equipment has been received in an undamaged condition, • All connections among instruments, computers, servers, and network are connected correctly, • Computers have correct CPU and memory, • Correct OS and software titles installed
OQ	<p>OPERATIONAL QUALIFICATION: <i>Documented verification that system components operate properly in accordance with manufacturer and customer requirements, government regulations, and industry standards. OQ establishes confidence that the components are capable of consistently operating within the established limits and tolerances required by the systems they support.</i></p>	<ul style="list-style-type: none"> • General system operations function correctly, • Audit trails capture correct information, • Electronic signatures are rendered correctly, • Security roles grant or restrict the correct permissions, • Users without credentials cannot log into system
PQ	<p>PERFORMANCE QUALIFICATION: <i>Documented verification that the total system performs as intended. Performance qualification establishes confidence that the system as a whole is capable of consistently performing within established limits.</i></p>	<ul style="list-style-type: none"> • Customer can go through typical laboratory workflow, • Customer's worst-case load scenario doesn't overwhelm server(s).

General guidelines regarding IQ, OQ, and PQ documentation

Aperio Spectrum Plus Operational Qualification Protocol

Customer Approval to Proceed with Protocol Execution

We have reviewed this protocol and prospectively approve it for qualification of the Spectrum Plus system described herein.

Printed Name and Title	Signature	Current Date
<i>Printed Name</i>		
<i>Title</i>		___/___/___
<i>Printed Name</i>		
<i>Title</i>		___/___/___

DSR Identification

<i>Aperio's System Model Number</i>	<i>Aperio's System ID</i>	<i>Recorded by: (Initials/Date)</i>

Customer Information

<i>Customer Name</i>	<i>Address of Installation Site</i>
<i>Customer Contact Person</i>	

Customer Approval of Completed Protocol

We have reviewed this protocol post-execution, and concur that operational qualification of the Spectrum Plus system described herein has been completed

Printed Name and Title	Signature	Current Date
<i>Printed Name</i>		
<i>Title</i>		___/___/___
<i>Printed Name</i>		
<i>Title</i>		___/___/___

FRM-0096	Aperio Spectrum Plus Operational Qualification Protocol	
Rev: 0	Aperio Technologies	Page 2 of 156
These Documents are the property of Aperio Technologies, Inc. and shall not be reproduced, distributed, disclosed or used for manufacture or sale of apparatus without the expressed written consent of Aperio Technologies, Inc.		

- Protocols must be prospectively prepared
- Protocols must be approved before being executed

General guidelines regarding IQ, OQ, and PQ documentation

- Protocols must be completed using good documentation practices, and documentation must occur at time of validation execution

Aperio Spectrum Plus Operational Qualification Protocol

7.3 Good Documentation Practices


This protocol shall be completed following Good Documentation Practices (GDP). Guidance below is based on Aperio procedure QSP-0023 (Ref 4.1.11). Consult the Customer to determine if the Customer has any additional documentation requirements that further govern documentation practices; these could further restrict the documentation instructions given below. In the event that Customer's instructions contradict those stated below, follow Customer's procedures.

7.3.1 Record all entries in indelible black or blue ink. Do not use ink of any other color, erasable ink, water-soluble ink, or pencil.


7.3.2 If recorded information must be changed, draw a single line through the incorrect recording, such that the original information is not obliterated. Record the new information as close as feasible to the original entry, and record your initials and the current date. If you think a reader may have any questions about the reason for the change, summarize the reason and record it next to the change.

Examples:


INCORRECT: Data is written over.

	List and sum the amperage of all other items drawing electrical current from this circuit	Sum recorded at right	6.5 AMPS
---	---	-----------------------	----------

INCORRECT: Data is obscured.

	List and sum the amperage of all other items drawing electrical current from this circuit	Sum recorded at right	5.5 AMPS 6.0
--	---	-----------------------	----------------------------

Correct documentation practices

	List and sum the amperage of all other items drawing electrical current from this circuit	Sum recorded at right	5.5 AMPS 6.0 ECS 01/01/08 Added incorrectly
---	---	-----------------------	---

7.3.3 When recording dates, use the format of MM/DD/YY, unless the Customer has another standard.

7.3.4 Do not leave any blank spaces in the protocol. If any fields in the protocol are purposefully left blank, indicate that they are not applicable by drawing a diagonal line through the field(s), writing "N/A," and recording your initials and the current date.

General guidelines regarding IQ, OQ, and PQ documentation

Step	Instructions	Expected Result	Actual Result	Initials/Date
29.	In the New Password and Retype New Password fields, enter a new password that contains at least one non-alphanumeric character and is at least 8 characters in length. Click Save .	Login completes	<input type="checkbox"/> As specified <input type="checkbox"/> Other (explain)	Pass Fail
30.	Click Log off	Spectrum displays Login Required screen	<input type="checkbox"/> As specified <input type="checkbox"/> Other (explain)	Pass Fail

Comments

Section results recorded by: _____ Date: _____
 Section results reviewed by: _____ Date: _____

- Results must be independently reviewed:
 - On a per-section basis
 - At conclusion of execution

Customer Approval of Completed Protocol


We have reviewed this protocol post-execution, and concur that operational qualification of the Spectrum Plus system described herein has been completed

Printed Name and Title	Signature	Current Date
<i>Printed Name</i>		____/____/____
<i>Title</i>		
<i>Printed Name</i>		____/____/____
<i>Title</i>		

Agenda

- General Considerations
 - **General Laboratory Practices**
 - history
 - agencies
 - guidelines
 - observations
 - **Comprehensive Solution to Regulatory Guidelines**
 - instrument qualification
 - **software validation**
- Advantages for Compliance


Data Repository/Viewer Validation



Processing Architecture

Technical challenges around Entire Slide Processing

- **Processing of huge Digital Images**
e.g. 15mm x 15mm specimen scanned at 20X yields an image with 30,000 x 30,000 pixels or 2.7 GB of data (~60MB compressed JPEG2000 file)
- **Server-Side Processing**
analysis runs server-side
- **Distributed Computing**



Securing the Digital Slide



- Raw data image files are stored in a secured folder in the DSR
- A checksum is calculated for each image file
- If corruption or alteration of the image occurs, Spectrum Plus will display an error message and will not open the slide

Managing Images and Data

Digital Pathology Information System

- Server storage, Web based access
- Customizable organization of images and information
- Supports many image types
- Software Integration
- Access Security
- Archiving & Retrieval
- User roles and permissions
- Remote viewing, collaborations, virtual review boards, etc
- Slide metadata
 - LIS Interface
 - Barcodes








Image Analysis


You use Stains to highlight Features
 e.g. - Cell Compartments
 (nuclei, membrane, cytoplasm)
 - Protein expressions





Where and how much staining is there?

Where and how many objects are there?
 e.g. - Tumor Cells

How much staining is there on different objects?
 e.g. - Cell Compartments





Multi-Color Images with RGB Color Pixels (Picture Elements)
 Pixels of different Colors (Stains)



Classify Pixels Measure Pixel by Color (Stain)(Stain) Intensities

Identify Objects by Color, Shape, & Size

Measure Object (Stain) Intensities



Validation: Comprehensive, Detailed, Documented

Aperio Data System Operational Qualification, SW Release 11	
Table of Contents	
Aperio Data System Operational Qualification, SW Release 11	
1	OPERATIONAL QUALIFICATION OVERVIEW
2	SCOPE
3	RESPONSIBILITIES
4	REFERENCES
5	DEFINITIONS AND ABBREVIATIONS
6	INSTRUCTIONS FOR PROTOCOL EXECUTION
6.1	Customer Notification
6.2	Protocol Pre-Approval
6.3	Good Documentation Practices
6.4	Protocol Execution
6.5	Per-Section Review
6.6	Protocol Deviations
6.7	Customer Approval of Completed Protocol
7	IDENTIFICATION OF EXECUTORS AND REVIEWERS
8	TEST MATERIALS AND PREPARATION FOR TESTING
8.1	Materials
8.2	Support Software
8.3	Protocol Execution Notes
9	GENERAL FUNCTIONALITY OF SPECTRUM PLUS
9.1	System Configuration
9.2	Spectrum Digital Slides, Specimens, and Regions
9.2.1	Digital Slides: Adding, Creating, Deleting, and Modifying
9.2.2	Creating, Deleting, and Modifying Regions
9.2.3	Creating, Deleting, and Modifying Regions
9.2.4	Creating, Deleting, and Modifying Regions
9.2.5	Additional Spectrum UI Verification
9.3	Editing Spectrum Plus Data Tables and Regions
9.3.1	Editing Data Tables
9.3.2	Population of Data Fields for Searches
9.3.3	Searches: Initiating, Saving, and Deleting
9.3.4	Search Accuracy
9.3.5	Return of Data Tables to Original State
10	GENERAL FUNCTIONALITY OF IMAGESCOPE SOFTWARE
10.1	Slide Viewing
10.1.1	Opening Slides and Viewing Related Information
10.1.2	Methods for Viewing Images
10.1.3	Image View Files (.SIS file)
10.1.4	Opening and Viewing Fluorescence Stained Images
10.2	Image Resolution
10.3	Image Adjustment and Management
10.3.1	Generalized Image Adjustment in ImageScope
10.3.2	Integrated Color Management (ICM)
10.3.3	Image Quality Enhancements ("IQ")
10.4	Annotations
10.4.1	Drawing of Shapes and Annotation Audit Trail
10.4.2	Accuracy of ImageScope Measurements
10.4.3	Annotation Layers
10.5	Linking
10.6	Tracking
10.7	Snapshots
10.8	Extraction of Regions
10.8.1	Saving Extracted Regions
10.8.2	Viewing Extracted Regions with ImageScope
10.8.3	Viewing Extracted Regions with External Viewing Software
10.9	Image Rotation
10.10	Smart Synchronization (SmartSync)
11	LOGICAL SECURITY
11.1	Account and Password Controls
11.1.1	Configuration of Account and Password Controls
11.1.2	Challenge of Account and Password Controls
11.1.3	Logging or Reporting Invalid Login Notifications
11.1.4	LDAP Account and Password Controls
11.2	Record Protection
11.2.1	Restrictions on Unauthorized Access to System and Files
11.2.2	Detection of Alteration to Image File
11.3	Roles
11.3.1	Assignment and Creation of Roles
11.3.2	Roles: Command Permissions
11.3.3	Roles: Data Table Field Permissions
11.4	Data Groups
11.5	Inactivity Timeouts
11.6	Communication Security
11.6.1	SSL
11.6.2	Access to Spectrum over the Internet
12	AUDIT TRAIL AND ELECTRONIC SIGNATURE
12.1	Configuration
12.1.1	General Audit Trail and Electronic Signature Configuration
12.1.2	Configuration of Audit Trail Reasons for Changes
12.1.3	Configuration of Hierarchy Data for Audit Trail and Electronic Signature Testing
12.2	Audit Trail Operations and Electronic Signature Executions
12.2.1	Controls on Signing
12.2.2	Audit Trail Operations: Capture of Reason for Change
12.2.3	Components, Controls, and Functionality for Electronic Signature
12.2.4	User Name Controls
12.2.5	Controls on Viewing and Generating Audit Trail Reports
12.2.6	Generation of Audit Trail Reports
12.2.7	Accuracy of Audit Trail Reports
12.2.8	Confirmation of Signer Identity
12.2.9	Permanence of Signatures and Audit Trail
12.2.10	Electronic Signature Functionality for LDAP-Authenticated Users
12.2.11	Removal of Status Vocabulary
13	WEBSCOPE
13.1	General WebScope Functionality
13.2	WebScope Annotations
14	ALGORITHMS
14.1	Configuration
14.2	Development of Macros from Algorithms; Functionality of Bright Field PPC Algorithm
14.3	Functionality of Customer-Purchased Bright Field Algorithms
14.4	Fluorescence Algorithms
14.5	Batch Analysis in Spectrum Plus
14.6	Audit Trail and Removal of Custom Macros from Database
15	DATA EXPORT
16	DIGITAL SLIDE CONFERENCING

Role-Based Security

- Administration control
- Logon required
- Password protected

- Role based permissions
 - Read/write
 - Data entry
 - Editing
 - Status changes
 - Data group
 - Study
 - Sharing
 - etc

User Details


Login: testuser123
Full Name: testuser123
Password: ●●●●●●●● (at least 5 characters)
Retype Password: ●●●●●●●●

User Permissions

Administrator: False


Data Group Name	Data Group Description	Access Level	Guest Access Level *
Dr. Colvins Data Group	Dr. Colvins Data Group	Read Only	None
Study 2	Study 2	No Access	None
Study 1	Study 1	No Access	None
HistoCore	HistoCore	No Access	None
Toxicology	Toxicology	No Access	None
Pharmacology	Pharmacology	No Access	None

Securing the Digital Slide



- Raw data image files are stored in a secured folder in the DSR
- A checksum is calculated for each image file
- If corruption or alteration of the image occurs, Spectrum Plus will display an error message and will not open the slide

Audit Trails



- Unlike the raw data image files, metadata can and will change as slides are annotated, reviewed, and processed.
- Metadata is maintained in a secure database where it cannot be manipulated from outside the Spectrum Plus application
- Changes to metadata are tracked in secure audit trail
- Critical data changes require reason for change and electronic signature

Audit Trails

Spectrum - Audit Trail - Windows Internet Explorer

http://localhost/AuditList.php?TableName=Project&AuditInserts&AuditUpdates&AuditDeletes&AuditStartDate=2011-06-29 00:00:00&AuditStopDate=2011-06-29 23:59:59&ids[]=4

File Edit View Favorites Tools Help

Spectrum -- Search -- Advanced Search User: Curtis, Role: Study_Pathologist aperio

Projects Specimens Digital Slides TMA Blocks Analysis SecondSlide Administrative Help Log off

Audit Trail

Legend:

- Insertions
- Updates
- Deletions
- Accesses

UserName: Curtis
 Date/Time: 2011-06-29 13:59
 Filters: Tables == Project: 4
 Audit Start Date: 2011-06-29 00:00:00
 Audit Stop Date: 2011-06-29 23:59:59
 Operations: Inserts, Updates, Deletes

Buttons: Generate PDF Report, Export To CSV File

Updated Project: 4

DateTime	User	Status
2009-03-06 14:08:05	hmeyer	In progress
2011-06-29 12:07:59	administrator	Completed

Updated Project: 4

DateTime	User	Status
2011-06-29 12:07:59	administrator	Completed
2011-06-29 12:13:41	administrator	In progress

Updated Project: 4 Authorized by Curtis Reason for Change: Changed Status

DateTime	User	Status
2011-06-29 12:13:41	administrator	In progress
2011-06-29 12:29:44	Curtis	Completed

Updated Digital Slide: 128 Reason for Change: New Results - repeat for kidney study

DateTime	User	Comment
2011-06-29 13:58:49	Curtis	Good results

Spectrum Plus + TMA Version 11.0.0.725 © Copyright 2006-2010 Aperio Technologies, Inc. All Rights Reserved.

Example protocol content: Spectrum Plus / DSR

- Equipment inventory (correct server and drive array received in good condition?)
- Site prepared?
- Server has correct components in place (CPU, RAM, network cards, etc.)?
- Connectivity complete?
- All software installed and configured correctly?

8.1.2 Multiple Storage Array (MSA)

Purpose	This section verifies that the HP MSA(s) and their components have been received in good condition and that they match the Customer's order			
Prerequisites and Comments	<ol style="list-style-type: none"> 1. Unpack the MSA(s). 2. Examine the MSA(s) and drives to determine if they have been received in good condition. 3. Slide out power supplies to determine part numbers and serial numbers. 4. Compare the received items to the customer's PO and/or PR in Attachment 2. 			

Component	Item Name or Part #	Serial #	Received in good condition?	Matches customer's order?
First MSA			Yes No	Yes No N/A
MSA's Left Side Power Supply			Yes No	Yes No N/A
MSA's Right Side Power Supply			Yes No	Yes No N/A
Additional MSA(s) (if applicable)			Yes No N/A	Yes No N/A
Additional MSA(s) Power Supplies (if applicable)			Yes No N/A	Yes No N/A

Hard drives within MSA(s):				
	Number of drives	GB each drive	Drives received in good condition?	Matches customer's order?
First MSA			Yes No	Yes No N/A
Additional MSA(s) (if applicable)			Yes No	Yes No N/A

Comments

Section results recorded by: _____ Date: _____

Section results reviewed by: _____ Date: _____

FRM-0094	Aperio Spectrum Plus / DSR Installation Qualification Protocol		Page 17 of 37
Rev. 2	Aperio Technologies		
<small>These Documents are the property of Aperio Technologies, Inc. and shall not be reproduced, distributed, disclosed or used for manufacture or sale of apparatus without the expressed written consent of Aperio Technologies, Inc.</small>			

Protocol content: Spectrum Plus / DSR

- General functionality of Spectrum Plus and the DSR
 - Creating, deleting, and modifying database objects, comments,
 - Attaching documents,
 - Database search and query functions.
- General functionality of ImageScope software
 - Annotation of slides; audit trail capture of changes,
 - Image adjustment, linking, tracking, extraction, etc,
 - Drawing of shapes and accuracy of linear and area measurement.
- Algorithms
- WebScope application
 - General functionality
 - Annotations in WebScope and audit trail
- Data export
- Digital slide conferencing
- Reporting
- Image file and database backup and restore
- GLP Archive utility*

Protocol content: Spectrum Plus / DSR

- General functionality of Spectrum Plus and the DSR
 - Creating, deleting, and modifying database objects, comments,
 - Attaching documents,
 - Database search and query functions.
- General functionality of ImageScope software
 - Annotation of slides; audit trail capture of changes,
 - Image adjustment, linking, tracking, extraction, etc,
 - **Drawing of shapes and accuracy of linear and area measurement.**
- Algorithms
- WebScope application
 - General functionality
 - Annotations in WebScope and audit trail
- Data export
- Digital slide conferencing
- Reporting
- Image file and database backup and restore
- GLP Archive utility*

Case Study: Spectrum Plus / DSR Validation of Linear and Area Measurement

Aperio ImageScope v10.0.16.1771 - [184.svs]

File Edit Image View Tools Window Help

Zoom: 8.0X

Annotations - Detailed View



Region	Length (um)	Area (um ²)	Text
1	956.3	57152	
2	780.3	48346	
3	177.8		

36447 x 32115 = 3.3GB, file = 64MB | 25605, 22460 : 3665 x 3760 | 26830, 22465 | prefetching / trackmap / progressive rendering | AAC PAN

How do we verify that these measurements of length and area are accurate?

Case Study: Spectrum Plus / DSR OQ Qualification of Linear and Area Measurement

Solution: The NIST-traceable Aperio Calibration Slide



Case Study: Spectrum Plus / DSR OQ

Validation of Linear and Area Measurement

Aperio Spectrum Plus Operational Qualification Protocol

Table for Recording Results for Accuracy of ImageScope Measurements

1. ImageScope will report perimeter and circumference in the Length (um) column of the Annotations Layer Regions table.
2. Calculation for Percent Difference is:
$$\frac{(100 * | \text{Value from Cal. Cert.} - \text{Value from ImageScope} |)}{\text{Value from Cal. Cert.}}$$


Object	Dimension	As Indicated on Cal. Cert.	As Calculated from Cal. Cert.	As Measured in ImageScope	Percent Difference
Large Rectangle (Nominally 15 mm wide x 6 mm high)	Width			From Step 11:	
	Height			From Step 11:	
Circle (Nominally 100 um diameter)	Width			From Step 13:	
	Circumference		Width x 3.1416 =	From Step 13:	
	Area		(Width/2) ² x 3.1416 =	From Step 13:	
Square (Nominally 10 um width)	Width			From Step 19:	
Square (Nominally 100 um height)	Height			From Step 16:	
	Perimeter		Height x 4 =	From Step 17:	
	Area		(Height) ² =	From Step 17:	

Comments

Presentation of IQ, OQ, PQ Results



SECTIONS AND CONTENTS OF THE SCANSCOPE AND SPECTRUM PLUS VALIDATION BINDER




TITLE PAGE AND INDEX

Contains:

- Description of the subject of the validation,
- Table of contents

Typical Size: 1 or 2 pages




VALIDATION FINAL REPORT

Contains:

- Summary of results from all IQ, OQ, PQ protocols
- Resolutions of any deviations found during protocol executions
- Signatures indicating that equipment is validated

Typical Size: 2 to 8 pages




VALIDATION PROTOCOLS: IQ, OQ, PQ

Contains:

- Completed IQ, OQ, and PQ protocols
- Signatures indicating that protocols were approved prior to execution
- Signatures indicating that work was reviewed after execution

Typical Size: 20 to 150 pages per protocol



EXHIBITS AND ATTACHMENTS

Contains:

- Printouts and other data generated during protocol execution
- Any supplemental testing required during protocol executions, for example, to further explore and resolve deviations
- Signatures indicating that data was reviewed after execution

Typical Size: 20 to 30 Exhibits of 1 to 10 pages each

Agenda

- General Considerations
 - **General Laboratory Practices**
 - history
 - agencies
 - guidelines
 - observations
 - **Comprehensive Solution to Regulatory Guidelines**
 - instrument qualification
 - software validation
- **Advantages for Compliance**

Aperio's GLP Compliance Module

Provide a customer-facing IQ/OQ/PQ package to accompany v11 release

- Why?
 - As noted, FDA and international regulatory agencies require qualification of such systems for all pre-clinical and clinical drug development
 - Ease customers' burden with compliance
 - Complements new regulatory-oriented functionality that comes with v11 release
 - Many other system/software manufacturers provide this type of service

Advantages of Using Vendor-Supplied Validation Services

- **Common in GxP laboratories and for GxP information systems.**
Examples include:
 - Networked chromatography systems and instrumentation (Agilent, Waters)
 - Laboratory Information Management Systems (SQL*LIMS, LabWare)
 - Electronic Document Management Systems (Livelink EDMS)
- **Advantages of using them:**
 - Vendor has comprehensive test case set that came from the development specifications of the hardware and software
 - Vendor's existing work saves time and cost of developing these on your own

Challenges in a vendor validation package

- Must be comprehensive but not excessive
- Must capture as much commonality as possible among various customer deployments
- Must assist with showing that other regulations are satisfied (e.g., ERES, HIPAA, CLIA, GCP, GMP, etc)

Digital Pathology Benefits for GLP

- **Archive** – breakage, photo-bleaching, oxidizing, loss, recovery, duplication, etc
- **Annotations** – Identify exact tumor mass/cell/organelle analyzed w/ scoring
- **Data Quality** – Objective data, blind studies, Image Analysis confirmation
- **Tracking** – All access, annotations, notes, image analyses, scoring, reports, etc are time/date stamped, tracked, electronically signed.....
- **Retrospective Analysis** – Raw image is always available for subsequent studies (review, new application, improved methods, new algorithms, etc)
- **Only practical method** to store fluorescent data
- **FDA submissions** – digital images easier to sort, search, transport, store, retrieve specific annotations and analyses

Acknowledgements

Thank you to the following Aperio team members for their assistance with this webinar:

- Craig Fenstermaker
- Mark Wrenn
- Kush Kapila
- Dr. Kate Lillard
- Jackie Gogue

And special thanks to our Compliance Consultant:

- Erik Staley, Valicom Inc.

References

- FDA regulations:
 - 21 CFR § 11 *Electronic Records; Electronic Signatures* (“Part 11”)
 - 21 CFR § 58, *Good Laboratory Practice for Nonclinical Laboratory Studies*
 - 21 CFR § 312, *Investigational New Drug Application* (Good Clinical Practices)
 - 21 CFR § 820 *Quality System Regulation* (medical devices and diagnostics)
 - 45 CFR § 160, 164, Standards for Privacy of Individually Identifiable Health Information (HIPAA)
- FDA Warning Letters: <http://www.fda.gov/foi/warning.htm>
- FDA and industry guidances
 - ICH, *Guidance for Industry, E6 Good Clinical Practice: Consolidated Guidance*
 - FDA Compliance Program Guidance Manual 7348.808, *Bioresearch Monitoring: Good Laboratory Practice (Nonclinical Laboratories)*
 - FDA, *Guidance for Industry, Computerized Systems Used in Clinical Investigations*
 - FDA, *Validation Documentation Inspection Guide*
 - FDA *Guidance for Industry, Part 11, Electronic Records; Electronic Signatures – Scope and Application*
 - FDA *Guidance for Industry, Quality Systems Approach to Pharmaceutical CGMP Regulations*
 - USP <1058>, *Analytical Instrument Qualification*
- EMEA regulations and guidances:
 - *OECD Principles of Good Laboratory Practice*, 1997
 - *The Application of the Principles of GLP to Computerised Systems*, 1995
 - Annex 11, “Computerised Systems,” 4/8/2008 draft, in *The Rules Governing Medicinal Products in the European Union*, Volume 4
- Japan MHLW regulations and guidances:
 - Ordinance No. 21, *Good Laboratory Practices* 3/26/97
 - PFSB Notification No. 0401022, “Use of Electromagnetic Records and Electronic Signatures for Approval of, or License for Drugs”
- Position papers:
 - *Toxicologic Pathology* 2007, 450-455, D. Tuomari et al, “Society of Toxicologic Pathology Position Paper on Pathology Image Data: Compliance with 21 CFR Parts 58 and 11

Questions

