

# ScanScope® XT / XT2 System

Experience ultra-fast, high-capacity scanning with our powerful 120-slide capacity ScanScope® XT scanner for 1 x 3" slides, or our 60-slide capacity XT scanner for 2 x 3" slides. Intended for high-throughput environments such as busy hospitals, reference labs, research labs, CROs, and pharma organizations, the ScanScope® XT and XT2 systems consistently deliver volumes of digital slides with superior image quality.

Aperio's high-capacity ScanScope® XT and XT2 systems deliver a high volume of digital slides with superior image quality.

After scanning your slides, you can view, manage, and analyze your data with our Spectrum™ digital pathology information management software.

## CREATE, VIEW, EDIT


- Rapidly load and scan with one-button automation
- Create digital slides in multiple formats (TIFF, composite web slide)
- View and edit digital slides with our user-friendly and freely downloadable ImageScope™ viewing software
- Connect multiple, remote parties through digital slide conferences

## MANAGE

- Manage digital slides, search data tables, and autopopulate fields with barcode data with Spectrum™, our web-based digital pathology information management software

## UPGRADE OPTIONS

- Spectrum™ Plus software for securely managing slides in multi-user, multi-site environments
- Optional 1-TB server highly recommended for optimal performance



# Specifications

## SLIDE LOADING

120 slide, automatic (1 x 3")

60 slide, automatic (2 x 3")

## SPEED (15 x 15mm)

<2 min/slide at 20x

## GLASS SLIDE DIMENSIONS

XT: 1 x 3" (25 x 75 mm) slides,  
with single 2 x 3" slide option

XT2: 2 x 3" (50 x 75 mm) slides,  
no 1 x 3" option

## OBJECTIVE

20x/0.75 Plan Apo

## SCANNING MAGNIFICATION

20x

40x (20x with 2x magnification changer)

## RESOLUTION

20x: 0.50 µm/pixel

40x: 0.25 µm/pixel

## IMAGE COMPRESSION

JPG2000, JPG

## IMAGE FILE FORMATS

TIFF (SVS), CWS, JP2

## SCANNING REGION (H x W)

23.4 x 55 mm

50 x 55 mm (large format)

## BARCODE SUPPORT

Interleaved DataMatrix

1-D 2-of-5

Code 39

Code 128

QR Code

## DIMENSIONS (W x D x H)

23.5 x 16 x 25.5" (597 x 406 x 648 mm)

## WEIGHT

105 lbs (47.6 kg)

## MONITOR

24" high-resolution LCD (1920 x 1200)

## POWER


AC 100-240, 50/60 Hz

## COMPLIANCE

CE, FCC

## WARRANTY

12 months



The ScanScope® XT features a 120-slide autoloader for high-throughput scanning, while the XT2 features a 60-slide autoloader.

For more information about our solutions for digital pathology, please visit [www.aperio.com](http://www.aperio.com) or call us at 866-478-4111.

**HEADQUARTERS**  
Aperio  
1360 Park Center Drive  
Vista, CA 92081

Main: 866.478.4111  
Fax: 760.539.1116

[info@aperio.com](mailto:info@aperio.com)  
[support@aperio.com](mailto:support@aperio.com)  
[www.aperio.com](http://www.aperio.com)

**EUROPEAN OPERATION**  
Aperio  
3 The Sanctuary  
Eden Office Park  
Ham Green  
Bristol BS20 0DD, UK

Main: +44 (0) 1275 375 123  
Fax: +44 (0) 1275 373 501

[europaero@aperio.com](mailto:europaero@aperio.com)  
[europesupport@aperio.com](mailto:europesupport@aperio.com)  
[www.aperio.com](http://www.aperio.com)



Make Pathology Extraordinary